

TITLE:	Policy and Procedure for the Acquisition and Use of Automated External Defibrillators (AEDs)	ROUTING Chief Operating Officer						
NUMBER:	BUL-4480.1	Educational Service Center: Superintendents Operation Administrators						
ISSUER:	Michelle King, Senior Deputy SuperintendentOperation CoordinSchool Operations, Office of the SuperintendentPrincipalsSchool PhysiciansSchool Physicians							
	René Gonzalez, Executive Director Student Health and Human Services	Nursing Coordinators School Nurses School Administrative						
DATE:	November 5, 2012	Assistants						
POLICY:	The Directors of Student Medical and Employee Health Services manage the Automated External Defibrillators (AEDs) Program of Los Angeles Unified School District (LAUSD). The policies and protocols provide uniformity and guidance in the administration and maintenance of the AED Program. Trained personnel will use an AED on persons who are unconscious, not breathing, or gasping, and not exhibiting signs of circulation. AEDs will be maintained on the premises of selected schools, and other locations in the Los Angeles Unified School District.							
MAJOR CHANGES:	This Bulletin replaces BUL-4480.0 of the same title dated December 19, 2008. The Annual Maintenance Inspection is the responsibility of the school administrator or designee. The Field Nursing Coordinator at each ESC is responsible for assisting with corrective actions on reported AED concerns and post-incident reviews. The Operations Director will ensure that schools are in compliance with Bulletin 4480.1.							
GUIDELINES :	: I. Responsibilities							
	A. Medical Director The Director of Student Medical Services will oversee all student related incidents, and the Director of Employee Health Services will oversee employees and all other adult related incidents. The Directors' responsibilities will be:							
	 To assist with the AED program and to act as consultant in matters regarding regulations and requirements for training and maintenance of AEDs in the AED program. 							
	2. To assist with policy development, procedures, and protocols for emergencies related to performing cardiopulmonary resuscitation (CPR) and use of an AED.							

3. To review recorded data from the AED Data Card or Data Port within five days following an AED incident involving a student or adult.



- B. AED Program Nursing Coordinator, representing the District Nursing Services, shall be responsible for the following:
 - 1. Coordinate, consult, and implement AED medical policies, protocols, CPR/AED training courses, and maintain training records.
 - 2. Coordinate the post-incident report and return the AED to readiness status.
- C. Educational Service Center (ESC) Nursing Field Coordinator
 - Assist with corrective action on reported AED problems in accordance with guidelines set forth by the manufacturer. The corrective actions will be documented on the AED Annual Maintenance Checklist (Attachment C) with the original kept at school site and copy sent to the respective ESC – Nursing Field Coordinator. The Annual Maintenance Checklist is to be completed annually and as needed.
 - 2. Participate in post-incidence reviews, debriefings, and quality assurance.
- D. Principal or Administrative Designee
 - 1. Ensure that all school administrators and staff annually receive a brochure (Attachment E) that describes the proper use of an AED and post the brochure next to every AED.
 - 2. Notify all employees annually of the location of all AED units on campus.
 - 3. Designate AED/CPR trained employees to be the AED Emergency Response Team, who shall be available to respond to an emergency that may involve the use of an AED. There should be a sufficient number of trained personnel to guarantee that at least one trained CPR/AED responder will be on site during school hours and school sponsored events.
 - 4. Confer with the AED Emergency Response Team members upon receipt of the AED, to determine:
 - a. A secure placement for the AED(s). Schools receiving their first AED should place it in an unlocked but alarmed cabinet located in the Main Office. All additional AEDs should be placed close to athletic events and practices, i.e., gymnasiums and fields.
 - b. Who will activate the site CPR/AED Emergency Response Team and how the Team will be notified.
 - c. Who will complete the CPR/AED Emergency Response Site Plan (Attachment A) and update annually. Mail a copy of Attachment A to



the AED Program Coordinator at District Nursing Services. The original copy is placed with the AED in the AED Program Protocol Book for three years.

- 5. Assign personnel to perform the Daily/Monthly Readiness Status Checklist (Attachment B) and perform the Annual Maintenance Checklist inspection (Attachment C). Mail a copy of Attachment C to the AED Program Coordinator each year in July. Original documents are placed with the AED in the AED Program and Protocol book for three years.
- 6. Schools hosting an athletic training or event should have an available AED. If a school has more than one AED unit and will be participating in an athletic event occurring at a non-LAUSD site, the designee should take an AED to that event. Please note all LAUSD high schools have AED(s).
- 7. Verify that CPR/AED cards are current for the Emergency Response Team. Cards from the American Heart Association are valid for two years, and cards from the American Red Cross are valid for one year.
- 8. Register for Heartsaver CPR/AED classes in the Learning Zone by clicking on the link: <u>http://lz.lausd.net/lz/index.jsp</u>. Click "Class Offerings" Go to "Keyword" and Type "CPR" Click Search.
- 9. It is recommended that schools consult with the AED Program Coordinator at District Nursing Services for the purchase of AEDs or prior to accepting any donated AEDs.
- E. AED Emergency Response Team
 - 1. Complete a basic CPR course that includes instruction on the use of the AED that complies with the regulations adopted by Emergency Medical Services authority (EMS) and standards of the American Heart Association (AHA) or the American Red Cross (ARC). The Heartsaver CPR/AED Course will satisfy this requirement.
 - 2. It is recommended that each member of the team have Hepatitis B vaccination. Hepatitis B vaccination is offered free of charge through LAUSD District Nursing Services. Call District Nursing for assistance regarding Hepatitis B vaccine.
- II. AED Equipment Maintenance and Supplies
 - A. AED equipment is the property of the school. The school will be responsible to reorder equipment supplies as noted below.



- 1. Automated External Defibrillator device case contents includes:
 - a. Automated External Defibrillator Device with battery and PC Data Card
 - b. Extra reserve battery in a sealed plastic covering
 - c. Two sets of adult defibrillator pads and one set of pediatric pads or pediatric key depending on the model
- 2. Fast Response Kit attached to AED Case includes:
 - a. Mouth-to-Mouth Protective Barrier/Mask
 - b. Scissors
 - c. One pair of disposable gloves
 - d. Disposable razor
 - e. Disposable towel
- B. Maintenance of AED Equipment

The AED located at the school or other site shall be maintained and regularly tested by the school principal or designee according to the operation and maintenance guidelines set forth by the manufacturer, and in accordance with any other applicable state and federal authority as noted below.

1. Daily AED Readiness Status Check

Trained Administrator/Designee will perform a Daily AED Readiness Status Check which includes the visual inspection of the AED to confirm the following:

- a. A flashing black hourglass means the FR2+ unit is ready to use. No Action is required.
- b. If a solid red "X" or flashing red "X" is seen in the Status Indicator Window, it should be reported for maintenance immediately. (For all units other than the FR2, check that the power indicator is present and the unit ready for use).
- c. For anything other than a flashing black hourglass or a solid/flashing red "X", remove and reinstall the battery to run the self-test. Follow manufacturer's maintenance task/response guidelines.
- d. Daily/Monthly Readiness Status Checklist (Attachment B) is completed daily.
- 2. Monthly Readiness Status Check Trained Administrator/Designee will perform Monthly AED readiness status check. The Monthly AED Readiness Status Check should confirm the following:
 - a. Defibrillator Pads: confirm that two sets of adult pads and one set Pediatric Pads are present and that they have not expired as indicated on the outside of the package.



- b. Spare Battery: present and has not expired as indicated on the battery "Install Before Date."
- c. Fast Response Kit: attached to the AED case and includes all the supplies listed under Equipment and Supplies.
- d. The Daily/Monthly Readiness Status Checklist (Attachment B) is completed.
- Annual Report of AED Maintenance Checklist Perform Maintenance Checklist inspection annually (Attachment C). Mail a copy of Attachment C to the Program Coordinator at District Nursing Services each year in July. Original documents are placed with the AED in the AED Program and Protocol book for three years.
- III. After Use of AED
 - A. The Principal or Administrative Designee will notify the District Nursing Services: AED Program Coordinator.
 - 1. The AED will be inspected and put back into readiness by the AED Program Coordinator/Designee from District Nursing Services and documented on the AED Program Daily/Monthly Readiness Status Checklist (Attachment B).
 - 2. <u>No AED unit should be given to any other agency until the unit is</u> inspected and data reviewed by the LAUSD Medical Director.
 - B. Complete the <u>Incident Report</u> in the <u>Incident System Tracking</u> <u>Accountability Report (ISTAR)</u> and the Confidential Report of Automated External Defibrillator (AED) Incident Response (Attachment D) for all events involving students, employees or visitors. Send both originals to the Director of the Office of Risk Management and Insurance Services and a copy to the Office of Environmental Health and Safety. *No copies of these reports are maintained at the school. For questions, the Administrator may contact the Office of the General Counsel.*
 - C. School nurses who respond to an incident that involves the use of an AED on a student will record their nursing assessment on the student's Health Record Card and/or in the Welligent Health Manager.
 - D. The Principal or Administrative Designee shall be responsible for organizing the post incident review to evaluate the effectiveness of the site AED response.
 - E. Order AED supplies, see Attachment F: The LifeTrends Group Phillips HeartStart Defibrillators.



AUTHORITY:	California Health and Safety Code §1797.196 California Civil Code § 1714.21 California AB 2083 California Code of Regulations Title 22. Social Security Division 9. Prehospital Emergency Medical Services Chapter 1.8
RELATED	Attachment A – CPR/AED Emergency Response Site Plan
RESOURCES:	Attachment B – Daily/Monthly Readiness Status Checklist
	Attachment C – AED Maintenance Checklist
	Attachment D – Confidential Report of AED Incident Response
	Attachment E – AED Information Pamphlet
	Attachment F – LifeTrends Group (LAUSD approved vendor). Order new equipment and supplies: (877) 443-2994.
	AED equipment repair: Philips Tech Support – (800) 263-3342
	American Heart Association
	American Red Cross
	Emergency Medical Services Authority
	OEHS – Safe School Plan, Volume 2 – Emergency Procedures
ASSISTANCE:	If you have any questions, please contact District Nursing Services, Special Programs/CPR/AED, at (213) 202-7580.



AUTOMATED EXTERNAL D	IFIED SCHOOL DISTRIC EFIBRILLATOR (AED) P CY RESPONSE SITE PL	ROGRAM
School Nurse:	Date:	
Educational Service Center: Loc Code:	School:	
School Address:		
City:	Zip	Code:
Phone () Extensio	n	
Exact Location of the AED		
CPR/AED Emergency Response Team Members	Employee No. CF	PR/AED Card Expiration Date

1. How will the first responder activate the Emergency Medical Service (EMS)?

- 2. How will the first responder notify the CPR/AED Emergency Response Team Members? Indicate by phone, radio, bell, intercom (i.e., "Code blue in room 20")
- 3. Which Emergency Response Team members have keys & access to the AED and will bring the AED to the emergency site?

The first responder will be responsible for documentation of the emergency.

- ✓ Please note Attachment A CPR/AED Emergency Response Site Plan must be reviewed and updated annually in July by the Principal or Administrative Designee.
- Send a copy to District Nursing Services. Attn: AED Program Coordinator. The AED Program Coordinator will forward a copy to the Medical Director.)
- ✓ Maintain the original copy of the CPR/AED Response Site Plan at your school site for 3 years.

Attachment A



											Α	UTO								LLATOR (atus Che	(AED) PROGRAM cklist	Attachme	ent B	
Scł	nool: _								N	Afr/M	odel:	<u>Phili</u>	ps He	artst	art FR	<u>2</u> Se	rial	No. (SN):		Location of AED):		
<u>DA</u> 1. `	Visual	lly in : Al	spect	the A iness-	ED for-us	se inc	iate da licator nainter	s, incl nance	luding					<i>'</i> "	M	тІХ	N	т	F	 Two a. 1 b. 0 2. Sparawithout 3. First 	LY CHECK : sets of defibrillator pads nspect for package for tamperir Check: 'Use Before' date indicat e battery check. Visually chec out opening the package. : aid emergency care kit	ed on the package.	e spare battery	
	ļ		h	ľ			ļ.	T H					T H			· ·		H			Disposable razor & scissors One pair of disposable gloves			
																				с. [Disposable towel Mouth-to-mouth protective barri			
																				4. Cheo	cking Readiness Status of AE	D Device		
																					Inspect AED readiness by view flashing green light in the status		he	
																				b. I	If a <u>solid</u> RED X or <u>flashing</u> RE	X is seen in the status indication	tor	
																					window, it should be reported for identified problem should be		onthly	
																				Actio	on Notification Log. Indicate whe			
																				probl 6. Save	lem. e original in the Program Prot	ocol Book for 3 years.		
																						-		
																				Caution: Electrical shock hazard. Dangerous high voltage & currents are present. <u>Do r</u> open the Philips FR2 (AED) unit, remove cover, or attempt repair. The Philips FA2 (AED				
																				unit shoul	d be returned to an authorized s			
																				<u>(800) 263</u>	<u>-3342</u> .			
Init	ial		Drint	Nan	סר								Si	anat	uro									
														-							MONTHLY VERIFICATIO	IN OF AED READINES	<u>8</u>	
																'use	befo	pad ore' dat Ped	te 'i	Spare battery install before' date	Print name	Signature	Date checked MM/DD/YY	
DA	TE/TI	IME	I	DEN	TIFIE	D PF	ROBL	EM –	- CON	NTAC	T PE	RSO	N											
													_											
													_											
													_											
													_											
	/																							



Attachment C

Inspected by _____ Date _____

AUTOMATED EXTERNAL DEFIBRILLATOR (AED) PROGRAM ANNUAL MAINTENANCE CHECKLIST

ESC	School
School Administr	ator
School Nurse	

Criteria for AED Mainter	namaa Chaak	Identified problem/Commente
Criteria for AED Mainter		Identified problem/Comments
AED # Serial Number:		
Indicate location of AED:		
Indicate type of storage:	_	
Standard AED Cabinet Yes N	0	
Other		
AED Locator Signs visible? Yes No		
Battery installed? Yes No		
Condition of battery Good Low		
Spare battery available? Yes No Indicate "Install Before Date"	:	
Status of visual alarm: Gradient Glass visible		
 Red X 		
Accurate date and time? Yes No		
Exterior components and sockets condition	Good 🗌 Poor	
Data card inserted? Yes No		
FORMS		Identified problem/Comments
Attachment A: CPR/AED Emergency Respon	se Site Plan	
Completed & updated Yes No Attachment B: Daily Readiness Status Check	liet	
Completed, signed, & dated I Yes I No		
Attachment B: Monthly Readiness Status Che Completed, signed, & dated Yes No	ecklist	
		Identified problem/Comments
Two sets of Adult AED pads in sealed pkg	☐ Yes ☐ No	
Adult pad Exp Date: 1) Exp Date		
One set of Pediatric AED pads in sealed pkg		
Pediatric pad Exp Date:		
FAST RESPONDER KIT available?		
1. Mouth-to-mouth protective barrier/mask		
2. Disposable razor	Yes No	
3. Scissors	□ Yes □ No	
4. One pair disposable gloves	Yes No	
5. Disposable hand towel	Yes No	

Other concerns: _____



Attachment D

LOS ANGELES UNIFIED SCHOOL DISTRICT CONFIDENTIAL REPORT OF AUTOMATED EXTERNAL DEFRIBRILLATOR (AED) INCIDENT RESPONSE

This is a confidential report for transmission to and use by attorneys for the Los Angeles Unified School District.

INSTRUCTIONS

- 1. Complete this form <u>Confidential Report of AED Incident Response</u> within 24 hours
- 2. Send this form to Office of Risk Management & Insurance Services; address listed below
- 3. Also complete and submit Incident Report within the Incident System Tracking Accountability Report (ISTAR)
- 4. No copy of AED Incident Response Report shall be retained by the school, or given to anyone, including the student or parent.

From:					
Name of School	Education Svc Center	Location Code	School Phone Number		
Name of Reporting Supervisor/Administrator	E-mail address		Date		
I. VICTIM INFORMATION					
	Home Addr	ess			
Last Name First Name	City	S	State Zip code		
Sex()M ()F Age Date of	Birth				
(Check one) Student Grade	LAUSD Employee #		Visitor		
II. AED INCIDENT SUMMARY					
Date of AED Incident Response	Time of incidentar	n/pm			
Exact location of incident					
Name of witness(s)					
Name of trained rescuer(s) responding					
(, , , , , , , , , , , , , , , , , , ,					
Emergency response site plan activated? Yes	No				
Was 9-1-1 called? Yes No If yes		9-1-1			
Was CPR given before the AED arrived? Yes					
Description of Incident					
Name of person completing form	Employe	ee # Dat	te		
ADDRESS AND CONTACT INFORMATION					
Office of Risk	Management & Insurance Ser	vices			

333 South Beaudry Avenue, 28th Floor Los Angeles, California 90017 Phone: (213) 241-3139 Fax: (213) 241-8993



LOS ANGELES UNIFIED SCHOOL DISTRICT POLICY BULLETIN

LAUSD District Nursing Special Programs CPR/AED American Heart Association Training Center

- Coordinates and oversees the LAUSD CPR/AED programs
- Offers courses for LAUSD employees:
 - ► CPR/AED training and recertification
 - ► First Aid classes

Visit the Learning Zone for classes or contact (213) 202-7580 for assistance.

https://lz.lausd.net/lz/index.jsp

CPR / AED RESOURCES

American Heart Association American Red Cross Emergency Medical Services Authority (EMS) National Center for Early Defibrillation U.S. Food and Drug Administration (FDA) BUL-4480.1 Attachment E

Rev. October 2012





Automated External Defibrillator (AED)

District Nursing Special Programs CPR/AED 121 N. Beaudry Avenue Los Angeles, CA 90012

> Phone: (213) 202-7580 Fax: (213) 580-6557 <u>http://dns.lausd.net</u>

Attachment E (1of 2)

LOS ANGELES UNIFIED SCHOOL DISTRICT

DISTRICT NURSING SERVICES Dee Apodaca, Director

Automated External Defibrillator (AED)



District Nursing Special Programs CPR/AED

Serop Hakimian Nursing Field Coordinator AED Program Coordinator



Background

Cardiovascular disease is the leading cause of death in the United States for both men and women. This disease has claimed the lives of more than 950,000 each year.

Sudden cardiac arrest (SCA) is the major complication of cardiovascular disease. Approximately 450,000 people suffer yearly from this event.

SCA is treatable. Treatment of SCA is an immediate shock to the heart, which stops the fatal rhythm and allows a normal heart rhythm to resume.

If the shock is performed in less than 3 minutes from onset, there is a 90% chance of survival. This shock can be performed by lay people or first responders by using an Automated External Defibrillator (AED).

The AEDs are easy, safe, and effective when used properly. A properly maintained AED and an AED action plan with trained personnel will provide immediate life-saving procedures to victims of SCA.

LAUSD CPR/AED PROGRAM

Onsite Location of the AED Unit(s):

Activate the Site Emergency Response Plan

- Determine if the scene is safe.
- Assess
 - If unresponsive, activate the Emergency Medical System (EMS).
- CALL 911
- Summon the AED Unit to be brought to the victim immediately.
- If not breathing, initiate CPR by a trained responder.

Operating an AED

- 1. **POWER ON the AED** (voice prompts instruct the rescuer).
- 2. ATTACH electrode pads to the victim's bare chest.
- 3. "CLEAR" the victim and allow the AED to ANALYZE the heart rhythm.
- 4. If SHOCK is advised, "CLEAR" the victim and push the SHOCK button when prompted and resume CPR.

If no SHOCK is advised, leave pads on victim and begin CPR.

- 5. After 2 minutes of CPR, the AED will reanalyze the rhythm. Follow the prompts to continue treating victim.
- Assist EMS personnel as directed until they take complete charge of the victim. The victim will be transported by EMS personnel.
- 7. The AED is to remain at the school or administrative office site.
- 8. Report the incident to the AED Program Coordinator as soon as possible or if after hours, notify the Local District Operations Administrator.

District Nursing Special Programs CPR/AED 121 N. Beaudry Avenue Los Angeles, CA 90012

Phone: (213) 202-7580 Fax: (213) 580-6557 http://dns.lausd.net



Attachment F





AUTHORIZED MASTER DISTRIBUTOR LAUSD Approved Vendor #189898

INFORMATION SHEET ORDERING WITHIN THE LAUSD DISTRICT PHILIP'S AED & ACCESSORIES (Automatic External Defibrillator)
AED Models within the District: • HeartStart FR2+ (M3861A) • HeartStart FRx (861304) • HeartStart Onsite (M5066A)
All Accessories: • Batteries • Pads – Child or Adult • AED Signs • Alarmed Cabinets/Keys • Fast Response Kits • FRx Child Key • Carrying Cases Please call: The LifeTrends Group, TLTG Inc. 1-877- 443-2994 <i>lifetrends3@cox.net</i> for Pricing / Ordering

34145 Pacific Coast Hwy. #101 Dana Point, CA 92629